

Product Datasheet

L-Selectin (orb372516)

Description	Peptide to L-Selectin
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Cys-Gln-Lys-Leu-Asp-Lys-Ser-Phe-Ser-Met-Ile-Lys
Formula	$C_{62}H_{105}N_{16}O_{18}S_2$
Purity	> 95%
Source	Synthetic
MW	1426.75
CAS Number	126880-86-2
Expiration Date	12 months from date of receipt.