

Product Datasheet

Lymphocyte Activating Pentapeptide (orb372517)

Description	Peptide to Lymphocyte Activating Pentapeptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Leu-Pro-Pro-Ser-Arg
Formula	$C_{25}H_{44}N_8O_7$
Purity	> 95%
Source	Synthetic
MW	568.67
CAS Number	120484-65-3
Expiration Date	12 months from date of receipt.

