

Product Datasheet

Lys-(Des-Arg9, Leu8)-Bradykinin (orb372521)

Description	Peptide to Lys-(Des-Arg9, Leu8)-Bradykinin
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Lys-Arg-Pro-Pro-Gly-Phe-Ser-Pro-Leu
Formula	$C_{47}H_{75}N_{13}O_{11}$
Purity	> 95%
Source	Synthetic
MW	998.20
CAS Number	71800-37-8
Expiration Date	12 months from date of receipt.

