

Product Datasheet

Lytic Peptide SB-37 (orb372526)

Description	Peptide to Lytic Peptide SB-37
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Met-Pro-Lys-Trp-Lys-Val-Phe-Lys-Lys-Ile-Glu-Lys-Val-Gly-Arg-Asn-Ile-Arg-Asn-Gly-Ile-Val-Lys-Ala-Gly-Pro-Ala-Ile-Ala-Val-Leu-Gly-Glu-Ala-Lys-Ala-Leu-Gly
Formula	$C_{188}H_{320}N_{54}O_{45}S$
Purity	> 95%
Source	Synthetic
MW	4089.05
Expiration Date	12 months from date of receipt.