

Product Datasheet

Mage-1 Antigen (161-169) (orb372529)

Description	Peptide to Mage-1 Antigen (161-169)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Glu-Ala-Asp-Pro-Thr-Gly-His-Ser-Tyr
Formula	$C_{41}H_{57}N_{11}O_{17}$
Purity	> 95%
Source	Synthetic
MW	975.97
CAS Number	144449-86-5
Expiration Date	12 months from date of receipt.

