

Product Datasheet

MAGE-3 Antigen (167-176) (orb372530)

Description	Peptide to MAGE-3 Antigen (167-176)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Met-Glu-Val-Asp-Pro-Ile-Gly-His-Leu-Tyr
Formula	$C_{53}H_{80}N_{12}O_{16}S$
Purity	> 95%
Source	Synthetic
MW	1173.36
CAS Number	178243-46-4
Expiration Date	12 months from date of receipt.

