

Product Datasheet

MBP (90-106) phosphorylated (orb372542)

Description	Peptide to MBP (90-106) phosphorylated
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Phe-Phe-Lys-Asn-Ile-Val-Thr-Pro-Arg-pThr-Pro-Pro-Pro-Ser-Gln-Gly-Lys-NH ₂
Formula	C ₈₉ H ₁₄₁ N ₂₅ O ₂₂
Purity	> 95%
Source	Synthetic
MW	1913.27
Expiration Date	12 months from date of receipt.