

Product Datasheet

Melanocyte Protein PMEL 17 (185-193) (orb372548)

Description	Peptide to Melanocyte Protein PMEL 17 (185-193)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ile-Thr-Asp-Gln-Val-Pro-Phe-Ser-Val
Formula	$C_{46}H_{72}N_{10}O_{15}$
Purity	> 95%
Source	Synthetic
MW	1005.14
CAS Number	162558-10-3
Expiration Date	12 months from date of receipt.

