

Product Datasheet

Melanocyte Protein PMEL 17 (256-264) (orb372549)

Description	Peptide to Melanocyte Protein PMEL 17 (256-264)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Tyr-Leu-Glu-Pro-Gly-Pro-Val-Thr-Ala
Formula	$C_{44}H_{67}N_9O_{14}$
Purity	> 95%
Source	Synthetic
MW	946.08
CAS Number	156761-76-1
Expiration Date	12 months from date of receipt.

