

Product Datasheet

Melanocyte Protein PMEL 17 (44-59) (orb372550)

Description	Peptide to Melanocyte Protein PMEL 17 (44-59)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Trp-Asn-Arg-Gln-Leu-Tyr-Pro-Glu-Trp-Thr-Glu-Ala-Gln-Arg-Leu-Asp
Formula	$C_{95}H_{137}N_{27}O_{29}$
Purity	> 95%
Source	Synthetic
MW	2105.32
Expiration Date	12 months from date of receipt.