

Product Datasheet

Met-Lys-Bradykinin (orb372553)

Description	Peptide to Met-Lys-Bradykinin
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Met-Lys-Arg-Pro-Pro-Gly-Phe-Ser-Pro-Phe-Arg
Formula	$C_{61}H_{94}N_{18}O_{13}S$
Purity	> 95%
Source	Synthetic
MW	1319.61
CAS Number	550-19-6
Expiration Date	12 months from date of receipt.

