

Product Datasheet

mini-ANP (orb372556)

Description	Peptide to mini-ANP
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Met-Cys-His-Phe-Gly-Gly-Arg-Met-Asp-Arg-Ile-Ser-Cys-Tyr-Arg-NH ₂ (Disulfide bridge Cys2-Cys13)
Formula	C ₇₅ H ₁₁₇ N ₂₇ O ₁₉ S ₄
Purity	> 95%
Source	Synthetic
MW	1829.19
Expiration Date	12 months from date of receipt.

