




Product Datasheet

LAG3 antibody (orb37257)

Description	LAG3 is a rabbit polyclonal antibody obtained from KLH conjugated synthetic peptide derived from hum...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide between 103-132 amino acids from the Central region of human LAG3.
Target	LAG3
Form/Appearance	Supplied in PBS with 0.09% (W/V) sodium azide
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Purity	This antibody is purified through a protein A column, followed by peptide affinity purification.
MW	57449 Da
Uniprot ID	P18627
NCBI	NP_002277.4
Dilution Range	WB~~1:1000, IHC-P~~1:10~50, FC~~1:25
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 cell line ...



Immunohistochemical analysis of formalin...