

Product Datasheet

Murine CMV pp 89 (170-174) (orb372572)

Description	Peptide to Murine CMV pp 89 (170-174)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	His-Phe-Met-Pro-Thr
Formula	$C_{29}H_{41}N_7O_7S$
Purity	> 95%
Source	Synthetic
MW	631.76
Expiration Date	12 months from date of receipt.

