

Product Datasheet

NES Topoisomerase II alpha (1017-1028) (orb372584)

Description	Peptide to NES Topoisomerase II alpha (1017-1028)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Asp-Ile-Leu-Arg-Asp-Phe-Phe-Glu-Leu-Arg-Leu-Lys
Formula	$C_{73}H_{117}N_{19}O_{19}$
Purity	> 95%
Source	Synthetic
MW	1564.86
Expiration Date	12 months from date of receipt.