

Product Datasheet

Neurogranin (28-43) (orb372585)

Description	Peptide to Neurogranin (28-43)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ala-Ala-Lys-Ile-Gln-Ala-Ser-Phe-Arg-Gly-His-Met-Ala-Arg-Lys-Lys
Formula	$C_{78}H_{134}N_{28}O_{19}S$
Purity	> 95%
Source	Synthetic
MW	1800.18
CAS Number	146554-17-8
Expiration Date	12 months from date of receipt.

