

Product Datasheet

Neuromedin (U25) (orb372590)

Description	Peptide to Neuromedin (U25)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Phe-Lys-Val-Asp-Glu-Glu-Phe-Gln-Gly-Pro-Ile-Val-Ser-Gln-Asn-Arg-Arg-Tyr-Phe-Leu-Phe-Arg-Pro-Arg-Asn-NH ₂
Formula	C ₁₄₄ H ₂₁₇ N ₄₃ O ₃₇
Purity	> 95%
Source	Synthetic
MW	3142.60
CAS Number	98395-76-7
Expiration Date	12 months from date of receipt.

