

Product Datasheet

Neuromedin (U8) (orb372591)

Description	Peptide to Neuromedin (U8)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Tyr-Phe-Leu-Phe-Arg-Pro-Arg-Asn-NH ₂
Formula	C ₅₄ H ₇₈ N ₁₆ O ₁₀
Purity	> 95%
Source	Synthetic
MW	1111.32
CAS Number	98395-75-6
Expiration Date	12 months from date of receipt.

