

Product Datasheet

Neuromedin S (orb372592)

Description	Peptide to Neuromedin S
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ile-Leu-Gln-Arg-Gly-Ser-Gly-Thr-Ala-Ala-Val-Asp-Phe-Thr-Lys-Lys-Asp-His-Thr-Ala-Thr-Trp-Gly-Arg-Pro-Phe-Phe-Leu-Phe-Arg-Pro-Arg-Asn-NH ₂
Formula	C ₁₇₃ H ₂₆₅ N ₅₃ O ₄₄
Purity	> 95%
Source	Synthetic
MW	3791.37
Expiration Date	12 months from date of receipt.