

Product Datasheet

Neuron Specific Peptide (orb372594)

Description	Peptide to Neuron Specific Peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Asp-Val-Ser-Asp-Gly-Ser-Ala-Glu-Arg-Arg-Pro-Tyr-Thr-Arg-Met-Gly-Ser-Gly-Gly-Leu-Lys-Leu-His-Cys-Val-His-Pro-Ala-Asn-Cys-Pro-Gly-Gly-Leu-Met-Val-Thr
Formula	$C_{161}H_{262}N_{52}O_{51}S_4$
Purity	> 95%
Source	Synthetic
MW	3870.44
Expiration Date	12 months from date of receipt.