

Product Datasheet

Neuropeptide AF (orb372595)

Description	Peptide to Neuropeptide AF
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ala-Gly-Glu-Gly-Leu-Asn-Ser-Gln-Phe-Trp-Ser-Leu-Ala-Ala-Pro-Gln-Arg-Phe-NH ₂
Formula	C ₉₀ H ₁₃₂ N ₂₆ O ₂₅
Purity	> 95%
Source	Synthetic
MW	1978.21
CAS Number	192387-38-5
Expiration Date	12 months from date of receipt.

