

Product Datasheet

Neuropeptide EI-Gly-Arg-Arg-MCH (orb372597)

Description	Peptide to Neuropeptide EI-Gly-Arg-Arg-MCH
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Glu-Ile-Gly-Asp-Glu-Glu-Asn-Ser-Ala-Lys-Phe-Pro-Ile-Gly-Arg-Arg-Asp-Phe-Asp-Met-Leu-Arg-Cys-Met-Leu-Gly-Arg-Val-Tyr-Arg-Pro-Cys-Trp-Gln-Val (Disulfide bridge Cys23-Cys32)
Formula	$C_{182}H_{282}N_{54}O_{52}S_4$
Purity	> 95%
Source	Synthetic
MW	4186.86
Expiration Date	12 months from date of receipt.