

Product Datasheet

Neuropeptide W-30 (orb372601)

Description	Peptide to Neuropeptide W-30
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Trp-Tyr-Lys-His-Val-Ala-Ser-Pro-Arg-Tyr-His-Thr-Val-Gly-Arg-Ala-Ser-Gly-Leu-Leu-Met-Gly-Leu-Arg-Arg-Ser-Pro-Tyr-Leu-Trp
Formula	$C_{165}H_{249}N_{49}O_{38}S$
Purity	> 95%
Source	Synthetic
MW	3559.17
CAS Number	383415-90-5
Expiration Date	12 months from date of receipt.