

Product Datasheet

Neuropeptide Y (Free acid) (orb372608)

Description	Peptide to Neuropeptide Y (Free acid)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Tyr-Pro-Ser-Lys-Pro-Asp-Asn-Pro-Gly-Glu-Asp-Ala-Pro-Ala-Glu-Asp-Met-Ala-Arg-Tyr-Tyr-Ser-Ala-Leu-Arg-His-Tyr-Ile-Asn-Leu-Ile-Thr-Arg-Gln-Arg-Tyr
Formula	$C_{189}H_{284}N_{54}O_{58}S$
Purity	> 95%
Source	Synthetic
MW	4272.80
Expiration Date	12 months from date of receipt.