

Product Datasheet

Obestatin (orb372618)

Description	Peptide to Obestatin
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Phe-Asn-Ala-Pro-Phe-Asp-Val-Gly-Ile-Lys-Leu-Ser-Gly-Val-Gln-Tyr-Gln-Gln-His-Ser-Gln-Ala-Leu-NH ₂
Formula	C ₁₁₆ H ₁₇₆ N ₃₂ O ₃₃
Purity	> 95%
Source	Synthetic
MW	2546.89
CAS Number	1081110-72-6
Expiration Date	12 months from date of receipt.