

Product Datasheet

p17 Gag (20-28) (orb372633)

Description	Peptide to p17 Gag (20-28)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Arg-Leu-Arg-Pro-Gly-Gly-Lys-Lys-Lys
Formula	$C_{45}H_{86}N_{18}O_{10}$
Purity	> 95%
Source	Synthetic
MW	1039.30
Expiration Date	12 months from date of receipt.