

## Product Datasheet

### p34 cdc2-derived peptide (orb372634)

|                         |                                                                     |
|-------------------------|---------------------------------------------------------------------|
| <b>Description</b>      | Peptide to p34 cdc2-derived peptide                                 |
| <b>Form/Appearance</b>  | Lyophilized powder                                                  |
| <b>Storage</b>          | Shipped at 4°C. Store at -20°C for one year.                        |
| <b>Note</b>             | For research use only                                               |
| <b>Protein Sequence</b> | Lys-Val-Glu-Lys-Ile-Gly-Glu-Gly-Thr-Tyr-Gly-Val-Val-NH <sub>2</sub> |
| <b>Formula</b>          | C <sub>62</sub> H <sub>104</sub> N <sub>16</sub> O <sub>19</sub>    |
| <b>Purity</b>           | > 95%                                                               |
| <b>Source</b>           | Synthetic                                                           |
| <b>MW</b>               | 1377.62                                                             |
| <b>Expiration Date</b>  | 12 months from date of receipt.                                     |