

Product Datasheet

p3K truncated (Lys58, Lys60, Lys63) Ea (54-68) (orb372637)

Description	Peptide to p3K truncated (Lys58, Lys60, Lys63) Ea (54-68)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Phe-Glu-Ala-Lys-Gly-Ala-Lys-Ala- Asn-Lys-Ala-Val-Asp
Formula	$C_{59}H_{97}N_{17}O_{19}$
Purity	> 95%
Source	Synthetic
MW	1348.53
Expiration Date	12 months from date of receipt.