

Product Datasheet

p60c-src Substrate I (orb372640)

Description	Peptide to p60c-src Substrate I
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Tyr-Ile-Tyr-Gly-Ser-Phe-Lys
Formula	$C_{44}H_{60}N_8O_{11}$
Purity	> 95%
Source	Synthetic
MW	877.01
Expiration Date	12 months from date of receipt.

