

Product Datasheet

p60c-src Substrate II Phosphorylated (orb372641)

Description	Peptide to p60c-src Substrate II Phosphorylated
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ac-Ile-pTyr-Gly-Glu-Phe-NH ₂
Formula	C ₃₀ H ₄₁ N ₆ O ₁₀ P
Purity	> 95%
Source	Synthetic
MW	676.65
CAS Number	284660-72-6
Expiration Date	12 months from date of receipt.

