

Product Datasheet

PACAP-38 (28-38) (orb372648)

Description	Peptide to PACAP-38 (28-38)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Gly-Lys-Arg-Tyr-Lys-Gln-Arg-Val-Lys-Asn-Lys-NH ₂
Formula	C ₆₁ H ₁₁₀ N ₂₄ O ₁₄
Purity	> 95%
Source	Synthetic
MW	1403.71
CAS Number	160489-86-1
Expiration Date	12 months from date of receipt.

