

Product Datasheet

PACAP-38 (6-38) (orb372650)

Description	Peptide to PACAP-38 (6-38)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Phe-Thr-Asp-Ser-Tyr-Ser-Arg-Tyr-Arg-Lys-Gln-Met-Ala-Val-Lys-Lys-Tyr-Leu-Ala-Ala-Val-Leu-Gly-Lys-Arg-Tyr-Lys-Gln-Arg-Val-Lys-Asn-Lys-NH ₂
Formula	C ₁₈₂ H ₃₀₀ N ₅₆ O ₄₅ S
Purity	> 95%
Source	Synthetic
MW	4024.80
CAS Number	143748-18-9
Expiration Date	12 months from date of receipt.