

Product Datasheet

Pancreatic Polypeptide (31-36) (orb372655)

Description	Peptide to Pancreatic Polypeptide (31-36)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Leu-Thr-Arg-Pro-Arg-Tyr-NH ₂
Formula	C ₃₆ H ₆₁ N ₁₃ O ₈
Purity	> 95%
Source	Synthetic
MW	804.00
Expiration Date	12 months from date of receipt.

