

Product Datasheet

Pancreatic Polypeptide (orb372656)

Description	Peptide to Pancreatic Polypeptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ala-Pro-Leu-Glu-Pro-Glu-Tyr-Pro-Gly-Asp-Asn- Ala-Thr-Pro-Glu-Gln-Met-Ala-Gln-Tyr-Ala-Ala- Glu-Leu-Arg-Arg-Tyr-Ile-Asn-Met-Leu-Thr-Arg- Pro-Arg-Tyr-NH ₂
Formula	C ₁₈₆ H ₂₈₇ N ₅₃ O ₅₆ S ₂
Purity	> 95%
Source	Synthetic
MW	4225.81
CAS Number	102082-88-2
Expiration Date	12 months from date of receipt.