

Product Datasheet

Pancreatic Polypeptide (Rana temporaria) (orb372657)

Description	Peptide to Pancreatic Polypeptide (Rana temporaria)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ala-Pro-Ser-Glu-Pro-His-His-Pro-Gly- Asp-Gln-Ala-Thr-Gln-Asp-Gln-Leu-Ala- Gln-Tyr-Tyr-Ser-Asp-Leu-Tyr-Gln-Tyr- Ile-Thr-Phe-Val-Thr-Arg-Pro-Arg-Phe- NH ₂
Formula	C ₁₉₂ H ₂₇₆ N ₅₂ O ₅₈
Purity	> 95%
Source	Synthetic
MW	4240.65
CAS Number	132187-74-7
Expiration Date	12 months from date of receipt.