

Product Datasheet

Pannexin-1 Fragment (4515) (orb372659)

Description	Peptide to Pannexin-1 Fragment (4515)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Glu-Lys-Asn-Ser-Arg-Gln-Arg-Leu-Leu-Asn-Pro-Ser
Formula	$C_{59}H_{104}N_{22}O_{20}$
Purity	> 95%
Source	Synthetic
MW	1441.62
Expiration Date	12 months from date of receipt.