

Product Datasheet

Pep 4A (orb372686)

Description	Peptide to Pep 4A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Gly-Ser-Val-Val-Ile-Val-Gly-Arg-Ile-Ile-Leu-Ser-Gly-Arg-NH ₂
Formula	C ₆₃ H ₁₁₇ N ₂₁ O ₁₆
Purity	> 95%
Source	Synthetic
MW	1424.77
Expiration Date	12 months from date of receipt.

