

Product Datasheet

PHI (orb372701)

Description	Peptide to PHI
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	His-Ala-Asp-Gly-Val-Phe-Thr-Ser-Asp-Phe-Ser-Arg-Leu-Leu-Gly-Gln-Leu-Ser-Ala-Lys-Lys-Tyr-Leu-Glu-Ser-Leu-Ile-NH ₂
Formula	C ₁₃₆ H ₂₁₆ N ₃₆ O ₄₀
Purity	> 95%
Source	Synthetic
MW	2995.40
CAS Number	80458-29-3
Expiration Date	12 months from date of receipt.

