

Product Datasheet

PHI (orb372702)

Description	Peptide to PHI
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	His-Ala-Asp-Gly-Val-Phe-Thr-Ser-Asp-Tyr-Ser-Arg-Leu-Leu-Gly-Gln-Ile-Ser-Ala-Lys-Lys-Tyr-Leu-Glu-Ser-Leu-Ile-NH2
Formula	$C_{136}H_{216}N_{36}O_{41}$
Purity	> 95%
Source	Synthetic
MW	3011.45
CAS Number	96849-38-6
Expiration Date	12 months from date of receipt.

