

## Product Datasheet

### PHM-27 (orb372703)

|                         |                                                                                                                             |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>      | Peptide to PHM-27                                                                                                           |
| <b>Form/Appearance</b>  | Lyophilized powder                                                                                                          |
| <b>Storage</b>          | Shipped at 4°C. Store at -20°C for one year.                                                                                |
| <b>Note</b>             | For research use only                                                                                                       |
| <b>Protein Sequence</b> | His-Ala-Asp-Gly-Val-Phe-Thr-Ser-Asp-Phe-Ser-Lys-Leu-Leu-Gly-Gln-Leu-Ser-Ala-Lys-Lys-Tyr-Leu-Glu-Ser-Leu-Met-NH <sub>2</sub> |
| <b>Formula</b>          | C <sub>135</sub> H <sub>214</sub> N <sub>34</sub> O <sub>40</sub> S                                                         |
| <b>Purity</b>           | > 95%                                                                                                                       |
| <b>Source</b>           | Synthetic                                                                                                                   |
| <b>MW</b>               | 2985.48                                                                                                                     |
| <b>CAS Number</b>       | 118025-43-7                                                                                                                 |
| <b>Expiration Date</b>  | 12 months from date of receipt.                                                                                             |

