

Product Datasheet

Phosphate Acceptor Peptide (orb372705)

Description	Peptide to Phosphate Acceptor Peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Arg-Arg-Lys-Ala-Ser-Gly-Pro-Pro-Val
Formula	$C_{41}H_{74}N_{16}O_{11}$
Purity	> 95%
Source	Synthetic
MW	967.15
Expiration Date	12 months from date of receipt.

