

Product Datasheet

Phosphorylase Kinase Beta-Subunit Fragment (420-436) (orb372707)

Description	Peptide to Phosphorylase Kinase Beta-Subunit Fragment (420-436)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Lys-Arg-Asn-Pro-Gly-Ser-Gln-Lys-Arg-Phe-Pro-Ser-Asn-Cys-Gly-Arg-Asp
Formula	$C_{79}H_{131}N_{31}O_{25}S$
Purity	> 95%
Source	Synthetic
MW	1946.20
CAS Number	150829-21-3
Expiration Date	12 months from date of receipt.

