

Product Datasheet

Phosphorylated Protein Kinase C Substrate 1 (orb372708)

Description	Peptide to Phosphorylated Protein Kinase C Substrate 1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Lys-Arg-Pro-pSer-Gln-Arg-His-Gly-Ser-Lys-Tyr-NH ₂
Formula	C ₅₇ H ₉₆ N ₂₃ O ₁₈ P
Purity	> 95%
Source	Synthetic
MW	1422.54
Expiration Date	12 months from date of receipt.

