

Product Datasheet

PKC Epsilon (orb372713)

Description	Peptide to PKC Epsilon
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Glu-Arg-Met-Arg-Pro-Arg-Lys-Arg-Gln-Gly-Ser-Val-Arg-Arg-Arg-Val
Formula	$C_{83}H_{155}N_{39}O_{21}S$
Purity	> 95%
Source	Synthetic
MW	2067.47
CAS Number	120253-69-2
Expiration Date	12 months from date of receipt.

