

Product Datasheet

PKI-tide (orb372714)

Description	Peptide to PKI-tide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ile-Ala-Ala-Gly-Arg-Thr-Gly-Arg-Arg-Gln-Ala-Ile-His-Asp-Ile-Leu-Val-Ala-Ala
Formula	$C_{85}H_{149}N_{31}O_{24}$
Purity	> 95%
Source	Synthetic
MW	1989.31
CAS Number	126370-52-3
Expiration Date	12 months from date of receipt.

