

Product Datasheet

Platelet Membrane Glycoprotein IIB Peptide (296-306) (orb372716)

Description	Peptide to Platelet Membrane Glycoprotein IIB Peptide (296-306)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Thr-Asp-Val-Asn-Gly-Asp-Gly-Arg-His-Asp-Leu
Formula	$C_{47}H_{75}N_{17}O_{20}$
Purity	> 95%
Source	Synthetic
MW	1198.22
Expiration Date	12 months from date of receipt.

