

Product Datasheet

Pneumadin (orb372720)

Description	Peptide to Pneumadin
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Tyr-Gly-Glu-Pro-Lys-Leu-Asp-Ala-Gly-Val-NH ₂
Formula	C ₄₇ H ₇₄ N ₁₂ O ₁₅
Purity	> 95%
Source	Synthetic
MW	1047.18
CAS Number	130918-90-0
Expiration Date	12 months from date of receipt.

