

Product Datasheet

Polylysine (orb372721)

Description	Peptide to Polylysine
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Lys-Lys-Lys-Lys-Lys-Lys-Lys-Lys-Lys
Formula	$C_{60}H_{122}N_{20}O_{11}$
Purity	> 95%
Source	Synthetic
MW	1299.77
CAS Number	25104-18-1
Expiration Date	12 months from date of receipt.