

Product Datasheet

pp60 c-SRC (521-533) Phosphorylated (orb372725)

Description	Peptide to pp60 c-SRC (521-533) Phosphorylated
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Thr-Ser-Thr-Glu-Pro-Gln-pTyr-Gln-Pro-Gly-Glu-Asn-Leu
Formula	$C_{62}H_{95}N_{16}O_{28}P$
Purity	> 95%
Source	Synthetic
MW	1543.53
CAS Number	149299-77-4
Expiration Date	12 months from date of receipt.

