

Product Datasheet

Preproenkephalin B (186-204) (orb372735)

Description	Peptide to Preproenkephalin B (186-204)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ser-Ser-Glu-Val-Ala-Gly-Glu-Gly-Asp-Gly-Asp-Ser-Met-Gly-His-Glu-Asp-Leu-Tyr
Formula	$C_{78}H_{115}N_{21}O_{36}S$
Purity	> 95%
Source	Synthetic
MW	1954.97
Expiration Date	12 months from date of receipt.

